

Microzone MRSA Detection Kit

Instructions For Use

MRSA-100

100 Reactions

Manual Date: October 2024

Version Number: 1.0

Location: Microzone,
Suite 3, Faraday House,
King William Street,
Stourbridge,
DY8 4HD

T: +44(0)1384 444585

E: info@microzone.co.uk

W: microzone.co.uk

Contents

Introduction.....	2
Kit Contents	2
Storage and Stability.....	2
Materials Required but Not Provided	2
Assay Principle	3
Technology	4
Assay Design	4
Protocol	5
Nucleic Acid Extraction/Purification.....	5
Reaction Setup.....	5
Thermal cycling.....	6
Interpretation of Results	6

Introduction

This kit is for the identification of Methicillin Resistant *Staphylococcus aureus* (MRSA) and other *Staphylococcus aureus* nucleic acids using Real-Time PCR (qPCR). It targets specific conserved regions of the of the bacteria's genome. It is intended for research use only by trained laboratory staff.

Kit Contents

Product	Product Code	Quantity
2X MegaMix Platinum Probe qPCR Mastermix	2MMPU-1	1 x 1 mL
MRSA Primer/Probe Mix	PP-MRSA-0.5	1 x 500 µL
Positive Control	PC-MRSA-0.1	1 x 100 µL
Just Water	5JWA-1	1 x 1 mL

Storage and Stability


Store the Microzone MRSA detection kit at -20°C until the expiry date displayed on the label.

Materials Required but Not Provided

- Disposable powder-free nitrile gloves
- Pipettes (adjustable) and sterile filtered pipette tips
- 0.5/1.5 mL nuclease free microcentrifuge tubes
- Benchtop centrifuge
- Vortex mixer
- Decontaminated bench
- Optically clear PCR tubes and caps or 96-well microplate suitable for use with chosen PCR instrument
- Plate sealer (if required)
- qPCR Thermal cycler

Assay Principle

The assay utilises qPCR to detect the presence of specific sequences of DNA found in target of interest. Specifically, the dual-labelled hydrolysis probe technology is employed as a method of detection as it allows the sensitive and specific detection of multiple targets within one reaction. The qPCR amplification consists of 4 phases throughout the run: initially it is a linear ground phase, where there is not sufficient amplified DNA for a change in fluorescence to be measured; next is an exponential phase, where detectable fluorescence is observed; followed by a Log-linear phase, where the efficiency of the reaction is measured; finally the reaction reaches the plateau phase, where the reaction slows down resulting in a reduced rate of increase in fluorescence as the reaction components deplete.



Linear plot of a qPCR curve.

Technology

Dual-labelled hydrolysis probes are the most widely employed method in quantitative PCR (qPCR) for detecting and quantifying specific DNA sequences. It utilises a fluorescently labelled probe that binds to a complementary target sequence within the PCR product. The probe is a short oligonucleotide, consisting of a reporter dye at the 5' end and a quencher at the 3' end. While the probe is intact, the quencher suppresses the fluorescence of the reporter. During the PCR amplification, the 5' to 3' exonuclease activity of Taq polymerase cleaves the probe's nucleotides, separating the reporter from the quencher, leading to an increase in fluorescence proportional to the amount of product amplified. The increase of fluorescence signal due to the dual-labelled hydrolysis probe technology allows the real-time monitoring of PCR amplification, providing highly specific and sensitive detection of nucleic acids.




Diagram showing how dual-labelled hydrolysis probes work.

Assay Design


This qPCR assay detects the *mecA* and *mecC* methicillin resistance genes of MRSA strains, and the species-specific *femA* and *nuc* genes of *Staphylococcus aureus*. It also incorporates an internal amplification control (IAC) which verifies that the conditions in each reaction are compatible with PCR amplification.

Protocol

The protocol consists of three steps: nucleic acid extraction/purification, qPCR and analysis. Nucleic acids are extracted/purified from samples suspected of being MRSA positive. The purified nucleic acids are then amplified using MegaMix Platinum qPCR Mastermix and MRSA Primer/Probe Mix. Analysis of the qPCR amplification allows the user to determine if the tested sample contained MRSA other *Staphylococcus aureus* nucleic acids.

Nucleic Acid Extraction/Purification

This kit can be used with tradition nucleic acid extraction methods that utilise magnetic beads or columns and alternatively, it can be used with a direct-to-PCR lysis buffer such as microLYSIS Plus from Microzone. microLYSIS Plus significantly shortens the sample preparation process, with only a short heat step to release and stabilise nucleic acids ready for PCR. More information about microLYSIS Plus can be found here: <https://microzone.co.uk/shop/extraction/microlysis-plus/>.



Reaction Setup

1. Thaw all reagents at room temperature.
2. Once thawed vortex all reagents and briefly centrifuge to remove any liquid from the cap.
3. Every reaction set up must allocate sufficient wells for the number of samples, a positive, and a negative control (samples + 2 = total wells needed). For example, 10 samples to test + 1 Positive control well + 1 negative control well = 12 reaction wells required.
4. Setup should be done on ice when possible. New sterile filter pipette tips are to be used between each pipetting.
5. Generate the following mastermix, use the number of reactions well determined in step 1, adding 5% excess (to allow for pipetting loss) to calculate volumes:

Component	Volume
MegaMix Platinum qPCR Mastermix with UNG	10 μ L
Primer and probe mix	5 μ L
Total volume of mastermix per sample	15 μL

6. Briefly vortex to mix.

7. Dispense 15 µL of mastermix into each reaction well.
8. Add 5 µL of sample to each well containing mastermix; for the positive control and negative control add 5 µL Positive Control and 5 µL Just Water, respectively.
9. Cap the tube with corresponding cap.
10. Briefly centrifuge to move all of the reaction to the base of the tube (some instruments this is not required as it is done during the PCR).
11. Transfer the tubes to the PCR instrument.

Thermal cycling

Run the following thermal cycling profile:

Temperature	Time	Cycles	Fluorescence collection	Step
95°C	2 minutes	1		Initial Denaturation
95°C	3 seconds	42		Denaturation
60°C	20 seconds		FAM, CY5, HEX	Annealing/Extension

Interpretation of Results

Target	Reporter Fluorophore	Sample	
		Cq Value	Interpretation
femA and nuc	FAM	<40	Positive
		≥40 or No amplification	Negative
mecA and mecC	Cy5	<40	Positive
		≥40 or no amplification	Negative
Internal Amplification Control (IAC)	HEX	<30	Positive
		≥30 or no amplification	Negative

	FAM	Cy5	HEX (IAC)	Reporting
<i>Positive Control</i>	+	+	+	Valid
<i>Negative Control</i>	-	-	+	Valid
<i>Case 1</i>	+	+	+/-	MRSA positive
<i>Case 2</i>	+	-	+/-	<i>Staphylococcus aureus</i> positive and MRSA negative
<i>Case 3</i>	-	-	+	<i>Staphylococcus aureus</i> and MRSA negative
<i>Case 4</i>	-	+	+	Inconclusive
<i>Case 5</i>	-	-	-	Invalid

Note: Any other combination of Positive or Negative Control result invalidates the run.